


Anaesthetic Gas Scavenging System



- GENTEC[®] Anaesthetic Gas Scavenging System is used for the safe removal of excess anaesthetic gas and anaesthetic vapors
- Complies with EN 737-2, EN 737-4 (TYPE 1)
- Color-coded for easy gas service identification according to ISO 32
- Flow: 25~50 l/min (Pressure drop: 0.15~0.29 psi)
- Air driving flow : 10~15 l/min (Pressure: 58 psi)
- 15 mm piping
- 100% tested and individually packaged